

**SIMPLE...
RELIABLE...
VALUE-PRICED**

PORTO-DENT MOBILE EVACUATION UNIT

Porto-Dent is perfect where central installation is impractical, and for practitioners with an extra office to furnish.

Compact

Measures only 12" x 15" x 26". It is no obstacle at chairside and stores inconspicuously. Shipping weight – 50 lbs.

Value Priced

Easily the least expensive of oral evacuation units. Excellent for limited budgets.

Ideal for Backup

A lifesaver when the central system breaks down. Doctor just rolls it out, plugs it in, and keeps operating.

Easy Maintenance

Non-corrosive collection canister cleans easily and is equipped with an automatic overflow shut-off. Motor comes with a one-year guarantee and can be replaced in minutes with a screwdriver.

Efficient

Adjustable motor speed can evacuate indefinitely at the rate of 75 cfm, or can be adjusted to the doctor's preference.

Special Uses

Needed especially by the dentist in a temporary location or when it is not feasible to install a central system. Perfect for such specialists as periodontists, endodontists, orthodontists, and hygienists.

Moves Easily

Rolls quietly on large carpet casters.

Finish

Ice White Matte or Walnut

Specifications

By-pass motor 7.5 Amps – 115V 60Hz.
Also available in 220/230V.




Includes two hoses and two mouthpieces. Shown with optional metal mouthpieces.

STERI-DENT

FROM
ALLIED
DIAGNOSTIC IMAGING RESOURCES, INC.
MAKERS OF
REDI-CHEM

STERI-DENT - Oral evacuation mouthpieces and accessories



400569 SV-55 — Dry Oral Funnel - made of sterile white polycarbonate

400549 OMP — Operative Mouthpiece, stainless 7/16" O.D.

400552 OSM — Surgical Mouthpiece - fits standard and saliva ejector adaptors

400579 OSS — Extra Sleeves for Saliva Ejector (Packed 12)

400577 OCB — Mouthpiece Cleaning Brushes - for cleaning oral mouthpieces

400596 OCB-L — Mouthpiece Cleaning Brush - for cleaning high volume mouth pieces and valves. Brush diameter 3/4", brush length 8"

101184 SV-57 — Mouthpiece Adaptor, Standard

400550 SV-58 — Mouthpiece Adaptor for Saliva Ejector (Fits into SV-57 adaptor)

(800) 346-DENT

All prices are FOB Factory, Leicester, NY. Prices and designs are subject to change without notice.



FROM
ALLIED
DIAGNOSTIC IMAGING RESOURCES, INC.
MAKERS OF
REDI-CHEM

Steri-Dent • 2364 Leicester Road, P.O. Box 175 • Leicester, NY 14481-0175

Toll Free: (800) 346-DENT in US • Phone: (585) 382-3223 • Fax: (585) 382-3736 • email: kwhite@cpac-fuller.com