

**SIMPLE...  
RELIABLE...  
VALUE-PRICED**

**STERI-DENT  
DRY HEAT  
STERILIZER**



UL  
LISTED

***Steri-Dent Dry Heat Sterilizes bring high quality, rugged performance, and total dependability to your practice.***



**Model # 200 – (two trays)**



**Model # 300 – (three trays)**

**The Best Method for Sharp Instruments**

Steri-Dent Dry Heat Sterilizers kill bacteria without moisture and will not dull sharp instruments. Plus, they can be used to dry wet packs or as a back-up sterilizer in your practice — the perfect complement to an autoclave.

**Safe and Simple to Operate**

Dry heat sterilization works without pressure, steam, or chemicals, which increases safety. Steri-Dent products feature a rapid-start heating element for fast operation.

**High-Quality**

Stainless steel construction means no rust — ever!

**Maintenance-Free**

Unlike wet heat methods, Steri-Dent Dry Heat Sterilizers require no routine cleaning.

**Convenient**

Utility trays can be sterilized with lids on and stored without contamination. Steri-Dent Sterilizers also have automatic timers with heat “hold” capability for continuous use in a busy practice.

**Value-Priced**

Low purchase cost and low operating cost — no distilled water to buy!



FROM  
**ALLIED**  
DIAGNOSTIC IMAGING RESOURCES, INC.  
A CINCINNATI COMPANY  
**MAKERS OF  
MEDICHEM**

# STERI-DENT - Sterilizer Specifications

Model	Trays	Tray Size		Outside dimensions			Inside (useful) space			Shipping Weight	Power Consumption
		W	H	W	H	D	W	H	D		
No. 200	2	10"	6 1/4"	15 1/2"	10 1/2"	9 1/2"	12 1/2"	5"	7 1/2"	17 lbs	500 Watts
No. 300	3	11 3/4"	7 1/4"	19"	14 1/2"	11"	15"	7 1/2"	9"	27 lbs	735 Watts

Operates on 115V AC electrical power. Inquire about other voltage.

## DRY HEAT SELF-SEALING BAGS

Use Dry Heat Sterilizing Bags as part of your infection control procedure for an added measure of safety. The azure strip turns brown when it has been exposed to high temperature. The paper used is specifically formulated for dry heat sterilizers.

2.5" x 1.5" x 10.5" (6.4cm x 3.8cm x 26.7cm)



## (800) 346-DENT

All prices are FOB Factory, Leicester, NY. Prices and designs are subject to change without notice.



FROM  
**ALLIED**  
DIAGNOSTIC IMAGING RESOURCES, INC.  
EST. 1982  
**MAKERS OF**  
**REDI-CHEM®**

Steri-Dent • 2364 Leicester Road, P.O. Box 175 • Leicester, NY 14481-0175

Toll Free: (800) 346-DENT in US • Phone: (585) 382-3223 • Fax: (585) 382-3736 • email: [kwhite@cpac-fuller.com](mailto:kwhite@cpac-fuller.com)