

PerfectView[®]

by CPAC Imaging



2364 Leicester Road
Leicester, New York 14481
Tel: 585-382-3223
Toll Free US and Canada: 1-800-333-XRAY (9729)
Fax: 585-382-9481
Web site: <http://www.thinkperfectview.com>
E-mail: sales@thinkperfectview.com



Orthopedic Illuminators

- 14 x 36" illuminator provides adequate viewing surface for longer films
- Baked enamel Hi-Reflective white finish
- Surface and recessed mounted available
- Automatic and standard switching options
- Optional film clips PV754

Deluxe Ultra Thinline Illuminators

Model Series 2000, 3000, and 4000

All PerfectView illuminators are constructed from 20-gauge cold rolled steel and finished inside and out with non-yellowing baked white enamel finish. Ensuring durability and providing the best possible viewing conditions has been our goal for over fifteen years.

Modern style and neutral color blends with any room decor. All model series are available for wall mount and recessed mounted applications.

The viewing area is one-piece white, translucent, shatterproof plastic 1/8" (3.2mm) thick. PerfectView illuminators feature easy to remove viewing panels that require no tools, for quick access to lamps.

Our state of the art light diffusing system provides the most uniform light distribution presently available in the "slim line" illuminator category.

- All illuminators except the Orthopedic model are UL and CUL approved
- Self-adjusting film grips
- Cluster switches as a standard feature
- All steel unweld construction
- Superior light diffusion system
- Baked enamel Hi-Reflective white finish
- Low leakage magnetic-NPF, UL and CSA approved ballasts
- 7500 hour daylight lamps at 750 lumens each
- Quick and friendly service



Deluxe Ultra Thinline Specifications

- All steel construction
- High luster, baked enamel finish
- Plexiglass® MC Acrylic Sheet viewing panel
- 110V power supply required

2000 Series—Two bulb, Thinline design deluxe illuminator. (1900 NITS)

3000 Series—Three bulb, Thinline design, deluxe illuminator (3000 NITS)

4000 Series—Four bulb, Thinline design, deluxe illuminator (3800 NITS)

Features:

- Units may be purchased individually or in single or double tier multiples, to a maximum of six viewing panels per tier
- Two, three, and four bulb units available
- Three way rocker switches on all automatic units
- Illuminators can be surface mounted or recessed into the wall

Options for 2000, 3000, and 4000 Series:

PV7000	Hospital Grade Plug
PV7100	Master Switch
PV7315-B	Metal 15 degree Angle Brackets for 2 Tier Units (Pair)
PV7330-A	Metal 30 degree Angle Brackets for 2 Tier Units (Pair)
PV7400	Handle Only
PV7401	Stands (2 per set)
PV7402	Combination (Handle and Stands)

Shipping Weights (in pounds)

	Wall Mount	Recessed		Wall Mount	Recessed
PV1436	33	77	PV3022	126	133
PV2001	18	46	PV3033	180	184
PV2002	33	72	PV3044	220	230
PV2003	44	97	PV3055	255	265
PV2004	108	115	PV3066	315	325
PV2005	120	130	PV4001	24	51
PV2006	150	157	PV4002	38	76
PV2011	76	86	PV4003	56	100
PV2022	123	127	PV4004	122	132
PV2033	165	169	PV4005	140	150
PV2044	210	215	PV4006	170	187
PV2055	225	234	PV4011	85	95
PV2066	288	320	PV4022	129	134
PV3001	21	48	PV4033	185	188
PV3002	36	74	PV4044	230	240
PV3003	53	100	PV4055	275	285
PV3004	118	127	PV4066	327	340
PV3005	130	140			
PV3006	160	165			
PV3011	80	90			

Viewing Area Options:



1 in 1




2 in 1



3 in 1



4 in 1




5 in 1




6 in 1




1 over 1



2 over 2




3 over 3



4 over 4



5 over 5



6 over 6


15° and 30° Bottom Angles




1 over 1




2 over 2




3 over 3



4 over 4



5 over 5



6 over 6

Specification Chart

For Wall or Desk Mounting Wall Mounted Dimensions

Viewing Area	W	H	D
14 x 36	15 7/8	39 1/4	3 1/2
14 x 17	15 7/8	20 1/4	3 1/2
2 in 1	29 13/16	20 1/4	3 1/2
3 in 1	43 7/8	20 1/4	3 1/2
4 in 1	57 15/16	20 1/4	3 1/2
5 in 1	71 15/16	20 1/4	3 1/2
6 in 1	86	20 1/4	3 1/2
1 over 1	15 7/8	40 11/16	3 1/2
2 over 2	29 13/16	40 11/16	3 1/2
3 over 3	43 7/8	40 11/16	3 1/2
4 over 4	57 15/16	40 11/16	3 1/2
5 over 5	71 15/16	40 11/16	3 1/2
6 over 6	86	40 11/16	3 1/2

For Recessed Mounting Actual Recessed Unit Size

Viewing Area	W	H	D
14 x 36	19	42 1/3	3 1/2
14 x 17	19	23 3/8	3 1/2
2 in 1	33	23 3/8	3 1/2
3 in 1	47	23 3/8	3 1/2
4 in 1	61 1/8	23 3/8	3 1/2
5 in 1	75 1/8	23 3/8	3 1/2
6 in 1	89 1/8	23 3/8	3 1/2
1 over 1	19	43 5/8	3 1/2
2 over 2	33	43 5/8	3 1/2
3 over 3	47	43 5/8	3 1/2
4 over 4	61 1/8	43 5/8	3 1/2
5 over 5	75 1/8	43 5/8	3 1/2
6 over 6	89 1/8	43 5/8	3 1/2

For Recessed Mounting Rough Opening for Recessed Mounting

Viewing Area	W	H	D
14 x 36	16 5/8	39 7/8	3 1/2
14 x 17	16 5/8	20 15/16	3 1/2
2 in 1	30 9/16	20 15/16	3 1/2
3 in 1	44 5/8	20 15/16	3 1/2
4 in 1	58 11/16	20 15/16	3 1/2
5 in 1	72 11/16	20 15/16	3 1/2
6 in 1	86 3/4	20 15/16	3 1/2
1 over 1	16 5/8	41 3/8	3 1/2
2 over 2	3 9/16	41 3/8	3 1/2
3 over 3	44 5/8	41 3/8	3 1/2
4 over 4	58 11/16	41 3/8	3 1/2
5 over 5	72 11/16	41 3/8	3 1/2
6 over 6	86 3/4	41 3/8	3 1/2

*Additional architectural information available. Please contact PerfectView@1-800-333-XRAY (9729)

Transfer Boards


Rigid polyethylene patient transfer board with built-in hand grips. For use in transferring patients from bed to stretcher or radiology table.

- Reduces patient discomfort
- Durable polyethylene design
- Easy to slide under patient
- Easy to clean
- May be sprayed with anti-static spray (F669)
- Contour hand grips (no sharp edges)
- Minimum distortion
- Reduces the risk of back injury

PV6101 (Standard) 23" W x 72" L x 3/10" D
 PV6102 (Deluxe) 23" W x 72" L x 3/10" D
 PV6103 MRT/CT 72" W x 19" L x 3/10" D



Deluxe




Standard

Corner Units


PerfectView will custom design corner modules to meet your space and design needs. Corner modules are available in white or wood grain laminate with 15-degree end panels for two tier units. Corner cabinets available are: single bank over bright spot; single bank over solid door for storage.

Call PerfectView for your custom quotation.
 1-800-333-XRAY (9729)



Transfer Board Wall Hanger

Solid hardwood with polyurethane coating
 (18" L x 5" H x 3/4" D)
 PV6104



About PerfectView

Perfectview illuminators and accessories are manufactured in the U.S. by CPAC Imaging, which offers a wide array of complementary products and services, including imaging chemicals and chemical mixers, processor cleaning and maintenance supplies, silver recovery equipment, cartridges, and silver refining services. For more information, including specification drawings for silver recovery units, please contact us at www.thinkperfectview.com. Shares of CPAC, Inc. stock are traded on the Nasdaq under the symbol: CPAK.

Statement of Warranty

All equipment is manufactured to exacting standards and has been tested and inspected for proper workmanship and performance before shipping.

Any parts which are defective will be repaired or replaced within a one year period after date of shipment, provided the equipment has been used according to the instruction manual and has not been abused or tampered with.

The company assumes no responsibility for damage in transit and the customer should present any claim for such damage to the carrier.

This warranty gives you specific legal rights. You may also have additional rights that vary from state to state.

Any unit to be repaired under warranty must be shipped, freight prepaid, or delivered to a facility authorized to render services provided hereunder. Returned unit must be either in its original package or a similar package affording an equal degree of protection. All units must have a Material Return Authorization code (MRA) visible on the returned item's product shipping label. MRAs can be obtained by calling 585-382-3223.

Acrylic viewing panels and lamps are not included under this warranty.

Repair, replacement, or other arrangements should be made through your PerfectView dealer representative.